# Parcan

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# Parcan

**3% stabilized solution of sodium hypochlorite for canal irrigation**

## Composition

3 % sodium hypochlorite, sodium chloride, sodium carbonate, solute qs ad 100 g.

## Indications

Canal irrigation for root canal preparation.

## Properties

The root canal preparation, essential procedure for the success of the obturation, needs the use of irrigant solutions combined with a chemo-mechanical reaming.

PARCAN, a 3% sodium hypochlorite stabilized solution meets the requirements of the expected effects of an irrigation solution, i.e. :

- chemical effects by dissolution of the organic matrix of the dentin, allowing to flare the canal easily.
- the bactericidal potency of sodium hypochlorite.

PARCAN is a stabilized and buffered sodium hypochlorite solution ready for use.

## Directions for use

Pump 10 ml of solution in the syringe of the PARCAN KIT ; then adjust the special needle provided with and carry out the canal irrigation.

## Warnings and precautions for use

So as to eliminate any risks of injecting the product beyond the apex it is strongly recommended :

- to check the apical closure (the use of the product for immature teeth or teeth with root resorption is contraindicated)
- to avoid using of the product under pressure.

Store away from heat and direct light.

## Presentation

Perfusion-type flask of 250 ml solution ready for use.

***For professional use only in the practice of dentistry.***

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