

## SAFETY DATA SHEET

### 1 - IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

**Identification of the substance or preparation:**

Name: DERMAFILM.

Product code: 0618W.

**Company/undertaking identification:**

Registered company name: SEPTODONT.

Address: 58 rue du Pont de Creteil.94100.Saint Maur des Fosses.France.

Telephone: +33(0)149767000. Fax:+33(0)148855401. Telex:.

msds@septodont.com

http://www.septodont.com

**Emergency telephone: +33(0)1.45.42.59.59.**

Association/Organisation: ORFILA .

### 2 - HAZARDS IDENTIFICATION

This product is classed: Highly flammable liquid.

This preparation is not classed as hazardous to health by directive 1999/45/EC.

### 3 - COMPOSITION/INFORMATION ON INGREDIENTS

Full text of risk phrases appearing in section 3: see section 16.

**Hazardous substances present on their own:**

(present in the preparation at a sufficient concentration to give it the toxicological characteristics it would have in a 100% pure state)

This preparation contains no hazardous substance in this category.

**Other substances representing a hazard:**

No known substance in this category present.

**Substances present at a concentration below the minimum danger threshold:**

INDEX	CAS	EC	Name	Symb.	R:	%
100808137	18472-51-0	242-354-0	CHLOREHEXIDINE DIGLUCONATE SOLUTION	Xi	36/37/38	2.5 <= x % < 10

**Other substances with occupational exposure limits:**

INDEX	CAS	EC	Name	Symb.	R:	%
603-002-00-5	64-17-5	200-578-6	ETHANOL	F	11	50 <= x % < 100

### 4 - FIRST AID MEASURES

As a general rule, in case of doubt or if symptoms persist, always call a doctor.

NEVER induce swallowing in an unconscious person.

### 5 - FIRE-FIGHTING MEASURES

Chemical powders, carbon dioxide and other extinguishing gas are suitable for small fires.

**Extinguishing media which must not be used for safety reasons:**

Sometimes, and due to the presence of organic material, a fire may produce a thick black smoke. Exposure to decomposition products may be hazardous to health.

Do not breathe in smoke.

Prevent the effluent of fire-fighting measures from entering drains or waterways.

**Special protective equipment for fire-fighters:**

Fire-fighting personnel are to be equipped with autonomous insulating breathing apparatus.

### 6 - ACCIDENTAL RELEASE MEASURES

**Personal precautions:**

Avoid inhaling the vapors.

Consult the safety measures listed under headings 7 and 8.

#### Environmental precautions:

- Contain and control the leaks or spills with non-combustible absorbent materials such as sand, earth, vermiculite, diatomaceous earth in drums for waste disposal.
- Prevent any material from entering drains or waterways.
- Use drums to dispose of waste recovered in accordance with applicable regulations (see heading 13).
- If the product contaminates waterways, rivers or drains, alert the relevant authorities in accordance with statutory procedures

#### Methods for cleaning up:

- Clean preferably with a detergent, do not use solvents.

## 7 - HANDLING AND STORAGE

The regulations relating to storage premises apply to workshops where the product is handled.

#### Handling:

- Prevent the formation of flammable or explosive concentrations in air and avoid vapor concentrations higher than the occupational exposure limits.

#### Fire prevention:

- Use the product in premises where there are no naked flames or other sources of ignition and have protected electrical equipment
- Keep packages tightly closed and away from sources of heat, sparks and naked flames.
- Do not use tools which may produce sparks. Do not smoke.
- Prevent access by unauthorised personnel.

#### Recommended equipment and procedures:

- For personal safety, see §8.
- Observe precautions stated on label and also industrial safety regulations
- Avoid inhaling vapors.

#### Prohibited equipment and procedures:

- Smoking, eating and drinking are prohibited in premises where the preparation is used
- Never open the packages under pressure

#### Storage:

- Keep the container tightly closed in a dry, well-ventilated place
- Keep away from all sources of ignition - do not smoke.
- Keep well away from all sources of ignition, heat and direct sunlight

## 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION

Use personal protection equipment as per Directive 89/686/EEC.

#### Exposure limit values per INRS ED 984:

France	VME-ppm:	VME-mg/m3:	VLE-ppm:	VLE-mg/m3:	Notes:	TMP N°:
64-17-5	1000	1900	5000	9500	-	84

#### Exposure limit values (2003-2006):

Germany/AGW	AGW:	AGW:	Faktor:	Bemerkungen:	
64-17-5	500 ml/m3	960 mg/m3	2(II)	DFG. Y	
UK/WELs	TWA:	STEL:	Ceiling:	Definition:	Criterion:
64-17-5	1000 ppm	-	-	-	-
UK/OES	TWA:	STEL:	Ceiling:	Definition:	Criterion:
64-17-5	1000 ppm	-	-	-	-
ACGIH/TLV	TWA:	STEL:	Ceiling:	Definition:	Criterion:
64-17-5	1000 ppm	-	-	-	-

#### Skin protection:

- For further information, see § 11 of S.D.S. - Toxicological information.

## 9 - PHYSICAL AND CHEMICAL PROPERTIES

#### General information:

- Physical state: fluid liquid

#### Important health, safety and environmental information:

- pH of the substance or preparation: not relevant.  
The pH is impossible to measure or its value is not relevant.
- Boiling point/boiling range: not specified.

Flash point interval:	Flash point: $\leq 21^{\circ}\text{C}$ .
vapour pressure:	not relevant.
Density:	not stated.
water solubility:	Dilutable.
<b>Other information:</b>	
melting point/melting range:	not specified.
Self-ignition temperature:	not specified.
Decomposition point/decomposition range :	not specified.

## 10 - STABILITY AND REACTIVITY

Any apparatus likely to produce a flame or to have a metallic surface at high temperature (burners, electric arcs, furnaces etc.) must not be allowed on the premises.

## 11 - TOXICOLOGICAL INFORMATION

No data is available regarding the preparation itself.

### Other data:

CAS 64-17-5 : IARC Groupe 1 (The agent is carcinogenic to humans).

## 12 - ECOLOGICAL INFORMATION

No ecological data on the product itself is available.

The product must not be allowed to run into drains or waterways.

## 13 - DISPOSAL CONSIDERATIONS

Do not pour into drains or waterways.

### Waste:

Recycle or dispose of waste in compliance with current legislation, preferably via a certified collector or company.

Do not contaminate the ground or water with waste, do not dispose of waste into the environment.

### Soiled packaging:

Empty container completely. Keep label(s) on container.

Give to a certified disposal contractor.

## 14 - TRANSPORT INFORMATION

Exempt from transport classification and labelling.

Transport product in compliance with provisions of the ADR for road, RID for rail, IMDG for sea and ICAO/IATA for air transport (ADR 2007 - IMDG 2006 - ICAO/IATA 2007).

## 15 - REGULATORY INFORMATION

This preparation was classified in compliance with the directive known as <All preparations> 1999/45/EC and its adaptations

In addition directive 2004/73/EC with the 29° adaptation of directive 67/548/EEC (Hazardous substances) have been taken into account.

This preparation is not classed as hazardous to health by directive 1999/45/EC.

### Preparation classification:



Highly flammable.

### Particular hazards associated with the preparation and safety recommendations:

R 11	Highly flammable.
S 16	Keep away from sources of ignition - No smoking.
S 9	Keep container in a well-ventilated place.

## 16 - OTHER INFORMATION

Since the user's working conditions are not known by us, the information supplied on this safety data sheet is based on our current level of knowledge and on national and community regulations.

The product must not be used for any purposes other than those specified under heading 1 without first obtaining written handling instructions.

It is at all times the responsibility of the user to take all necessary measures to comply with legal requirements and local regulations.

The information given on this safety data sheet must be regarded as a description of the safety requirements relating to our product and not a guarantee of its properties

**Full text of risk phrases appearing in section 3:**

R 11 Highly flammable.

R 36/37/38 Irritating to eyes, respiratory system and skin.

**Labelling for contents (REGULATION (EC) No 648/2004 - 907/2006):**

- disinfectants