

## SAFETY DATA SHEET

### 1 - IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND OF THE COMPANY/UNDERTAKING

#### Identification of the substance or preparation:

Name: TEMPOPHORE.

Product code: 252SF.

#### Company/undertaking identification:

Registered company name: SEPTODONT.

Address: 58 rue du Pont de Creteil.94100.Saint Maur des Fosses.France.

Telephone: +33(0)149767000. Fax: +33(0)148855401.

msds@septodont.com

http://www.septodont.com

#### Emergency telephone: +33 (0)1 45 42 59 59.

Association/Organisation: INRS / ORFILA <http://www.centres-antipoison.net>.

### 2 - HAZARDS IDENTIFICATION

This product is not classed as flammable. Refer to the recommendations regarding the other products present on the site

Possibility of corrosive effects.

Possibility of noxious effects with symptoms of slight poisoning by contact with the skin and swallowing.

Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

#### Preparation classification:



Corrosive



Dangerous for the environment

R 50/53

Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

R 21/22

Harmful in contact with skin and if swallowed.

R 34

Causes burns.

R 44

Risk of explosion if heated under confinement.

For professional use only.

#### Other data:

Packaging < 125 ml.

### 3 - COMPOSITION/INFORMATION ON INGREDIENTS

Full text of risk phrases appearing in section 3: see section 16.

#### Hazardous substances present on their own:

(present in the preparation at a sufficient concentration to give it the toxicological characteristics it would have in a 100% pure state)

INDEX	CAS	EC	Name	Symb.	R:	%
10808019	1314-13-2		ZINC OXIDE	N	50/53	25 <= x % < 50
10808027	89-83-8		THYMOL	C	34 22	10 <= x % < 25

#### Other substances representing a hazard:

INDEX	CAS	EC	Name	Symb.	R:	%
10808016	8021-39-4		CREOSOTE	T N	34 50 24/25	2.5 <= x % < 10
10808023	75-47-8		IODOFORM	Xn	20/21/22 36/37/38 44	10 <= x % < 25
10808026	244-350-4		CAMPHRE RACEMIQUE	Xn	11 36/37/38 22	2.5 <= x % < 10
10808043	2216-51-5		MENTHOL CRISTALLISE	Xi	38	2.5 <= x % < 10

#### Substances present at a concentration below the minimum danger threshold:

No known substance in this category present.

#### Other substances with occupational exposure limits:

No known substance in this category present.

### 4 - FIRST AID MEASURES

As a general rule, in case of doubt or if symptoms persist, always call a doctor.

NEVER induce swallowing in an unconscious person.

**In the event of splashes or contact with eyes:**

Wash thoroughly with soft, clean water for 15 minutes holding the eyelids open.  
If there is any redness, pain or visual impairment, consult an ophthalmologist.

**In the event of splashes or contact with skin:**

Remove all contaminated or stained clothing immediately. Do not use them again until they have been decontaminated.  
Wash the skin thoroughly with soap and water or a recognised cleaner.  
DO NOT use solvents or thinners.  
If the contaminated area is widespread and/or there is damage to the skin, a doctor must be consulted or the patient transferred to hospital

**In the event of swallowing:**

In the event of swallowing, if the quantity is small (no more than one mouthful), rinse the mouth with water, administer activated medical charcoal and consult a doctor.  
Consult a doctor and show him the label.  
If swallowed accidentally, do not allow to drink, do not induce vomiting and transfer to hospital immediately by ambulance. Show the label to the doctor.

---

**5 - FIRE-FIGHTING MEASURES**

Not relevant.

---

**6 - ACCIDENTAL RELEASE MEASURES**

**Personal precautions:**

Avoid any contact with the skin and eyes.  
Consult the safety measures listed under headings 7 and 8.

**Environmental precautions:**

Prevent any material from entering drains or waterways.  
Use drums to dispose of waste recovered in accordance with applicable regulations (see heading 13).  
If the product contaminates waterways, rivers or drains, alert the relevant authorities in accordance with statutory procedures

**Methods for cleaning up:**

Clean preferably with a detergent, do not use solvents.

---

**7 - HANDLING AND STORAGE**

The regulations relating to storage premises apply to workshops where the product is handled.

**Fire prevention:**

Prevent access by unauthorised personnel.

**Recommended equipment and procedures:**

For personal safety, see §8.  
Observe precautions stated on label and also industrial safety regulations  
Product must not come into contact with the skin and eyes  
Provide safety showers and eye washes in workshops where the product is handled constantly

**Prohibited equipment and procedures:**

Smoking, eating and drinking are prohibited in premises where the preparation is used  
Never open the packages under pressure

**Storage:**

Keep the container tightly closed in a dry, well-ventilated place  
Keep away from food and drink, including those for animals.  
Keep away from heat and in a cool place.

---

**8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

Use personal protection equipment as per Directive 89/686/EEC.

**Exposure limit values per INRS ED 984:**

France	VME-ppm:	VME-mg/m3:	VLE-ppm:	VLE-mg/m3:	Notes:	TMP N°:
75-47-8	0.6	10	-	-	-	-
1314-13-2	-	5	-	-	-	-

**Exposure limit values (2003-2006):**

UK/WELs	TWA:	STEL:	Ceiling:	Definition:	Criterion:
75-47-8	0.6 ppm	1 ppm	-	-	-

UK/OES	TWA:	STEL:	Ceiling:	Definition:	Criterion:
75-47-8	0.6 ppm	1 ppm	-	-	-
1314-13-2	5 mg/m3	10 mg/m3	-	-	-
ACGIH/TLV	TWA:	STEL:	Ceiling:	Definition:	Criterion:
75-47-8	0.6 ppm	-	-	-	-
1314-13-2	2 mg/m3	10 mg/m3	-	-	R

**Hand protection:**

Protective creams may be used for exposed skin, but they should not be applied after contact with the product.  
Where there is a risk of contact with the hands, suitable gloves must be used

**Eye and face protection:**

Avoid contact with the skin and eyes  
Provide personnel with gloves, face masks and safety goggles.  
Provide eye washes in workshops where the product is constantly handled

**Skin protection:**

Wear suitable protective clothing, in particular an apron and boots. These items must be kept in good condition and cleaned after use.  
For further information, see § 11 of S.D.S. - Toxicological information.

---

**9 - PHYSICAL AND CHEMICAL PROPERTIES**

**General information:**

Physical state: paste

**Important health, safety and environmental information:**

pH of the substance or preparation: not relevant.

The pH is impossible to measure or its value is not relevant.

Boiling point/boiling range: not specified.

Flash point interval: not relevant.

Vapour pressure: not relevant.

Density: not stated.

Water solubility: Insoluble.

**Other information:**

melting point/melting range: not specified.

Self-ignition temperature: not specified.

Decomposition point/decomposition range : not specified.

---

**10 - STABILITY AND REACTIVITY**

This preparation, not classed as explosive, may nevertheless display explosive properties in practise if it is heated in a sufficiently confined environment.

---

**11 - TOXICOLOGICAL INFORMATION**

No data is available regarding the preparation itself.

It is generally agreed that substances contained are likely to cause tissue damage in less than four hours after application to the healthy, unbroken skin of an animal.

---

**12 - ECOLOGICAL INFORMATION**

No ecological data on the product itself is available.

The product must not be allowed to run into drains or waterways.

**Ecotoxicity:**

Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

**In accordance with Directive 2006/8/EC:**

CAS	EC	
1314-13-2		ZINC OXIDE

---

**13 - DISPOSAL CONSIDERATIONS**

Do not pour into drains or waterways.

**Waste:**

Recycle or dispose of waste in compliance with current legislation, preferably via a certified collector or company.

Do not contaminate the ground or water with waste, do not dispose of waste into the environment.

**Soiled packaging:**

Empty container completely. Keep label(s) on container.  
 Give to a certified disposal contractor.

**14 - TRANSPORT INFORMATION**

Transport product in compliance with provisions of the ADR for road, RID for rail, IMDG for sea and ICAO/IATA for air transport (ADR 2009 - IMDG 2008 - ICAO/IATA 2009).

Classification:



Pollutants to the aquatic environment:



UN1759=CORROSIVE SOLID, N.O.S.  
 (thymol)

ADR/RID	Class	Code	Pack gr.	Label	Ident.	LQ	Provis.	EQ	Cat.	Tunnel
	8	C10	III	8	80	LQ24	274	E1	3	E

IMDG	Class	2°Label	Pack gr.	LQ	EMS	Provis.	EQ
	8	-	III	5 kg	F-A,S-B	223 274	E1

IATA	Class	2°Label	Pack gr.	Passager	Passager	Cargo	Cargo	note	EQ
	8	-	III	822	25 kg	823	100 kg	A3	-
	8	-	III	Y822	5 kg	-	-	A3	-

**15 - REGULATORY INFORMATION**

This preparation was classified in compliance with the directive known as <All preparations> 1999/45/EC and its adaptations  
 In addition directive 2004/73/EC with the 29° adaptation of directive 67/548/EEC (Hazardous substances) have been taken into account.  
 This product is not classed as flammable.  
 Packaging < 125 ml.

**Preparation classification:**



Corrosive



Dangerous for the environment

**Contains:**

CREOSOTE  
 IODOFORM  
 CAMPHRE RACEMIQUE  
 THYMOL

**Particular hazards associated with the preparation and safety recommendations:**

- R 50/53 Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
- R 21/22 Harmful in contact with skin and if swallowed.
- R 34 Causes burns.
- R 44 Risk of explosion if heated under confinement.
- S 26 For professional use only.
- S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- S 36/37/39 Wear suitable protective clothing, gloves and eye/face protection.
- S 45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
- S 61 Avoid release to the environment. Refer to special instructions/Safety data sheets.
- S 60 This material and its container must be disposed of as hazardous waste.

S 57 Use appropriate container to avoid environmental contamination.  
Product submitted to restrictions on the marketing and use: see Directive 76/769/EEC and its ensuing adaptations.

**Particular provisions:**

Product not to be sold to the general public. (Dir 76/769/EC) - For professional use only.

---

**16 - OTHER INFORMATION**

Since the user's working conditions are not known by us, the information supplied on this safety data sheet is based on our current level of knowledge and on national and community regulations.

The product must not be used for any purposes other than those specified under heading 1 without first obtaining written handling instructions.

It is at all times the responsibility of the user to take all necessary measures to comply with legal requirements and local regulations.

The information given on this safety data sheet must be regarded as a description of the safety requirements relating to our product and not a guarantee of its properties

**Full text of risk phrases appearing in section 3:**

R 11	Highly flammable.
R 20/21/22	Harmful by inhalation, in contact with skin and if swallowed.
R 22	Harmful if swallowed.
R 24/25	Toxic in contact with skin and if swallowed.
R 34	Causes burns.
R 36/37/38	Irritating to eyes, respiratory system and skin.
R 38	Irritating to skin.
R 44	Risk of explosion if heated under confinement.
R 50	Very toxic to aquatic organisms.
R 50/53	Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.