



**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

- **1.1 Product identifier**
- **Trade name: DETARTRINE Z**
- **Article number / formula: 138S**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the preparation : dental product**
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
SEPTODONT.  
58 rue du Pont de Creteil  
94100 Saint Maur des Fosses  
France  
Téléphone: +33(0)149767000  
Fax: +33(0)148855401  
msds@septodont.com  
http://www.septodont.com
- **1.4 Emergency telephone number:**  
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
Belgique - Tel : 32 070/245 245  
UK : +44 (0)-870-190-6777 (CARECHEM24)

**SECTION 2: Hazards identification**

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**  
The product is not classified according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **Additional information:**  
EUH208 Contains Menthone. May produce an allergic reaction.  
EUH210 Safety data sheet available on request.
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

**SECTION 3: Composition/information on ingredients**

- **3.2 Chemical characterisation: Mixtures**
- **Description:** Mixture: consisting of the following components.

(Contd. on page 2)

GB

## Safety data sheet

### According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 1)

**Dangerous components:**

CAS: 10101-52-7	Zirconium(IV) silicate ⚠ Acute Tox. 4, H332; Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335	10-25%
CAS: 64-17-5 EINECS: 200-578-6 Index number: 603-002-00-5 RTECS: KQ 6300000 Reg.nr.: 01-2119457610-43-xxxx	ethanol ⚠ Flam. Liq. 2, H225; ⚠ Eye Irrit. 2, H319	≤2.5%

· SVHC No

 · **Additional information:** For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **General information:** No special measures required.
- **After excessive inhalation:** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult a doctor.
- **4.2 Most important symptoms and effects, both acute and delayed :**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed :**  
No further relevant information available.

### SECTION 5: Firefighting measures

- **5.1 Suitable extinguishing agents:**  
CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

### SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :** Not required.
- **6.2 Environmental precautions:** Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.
- **6.4 Reference to other sections :**  
No dangerous substances are released.  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protective equipment.  
See Section 13 for disposal information.

### SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :** No special measures required.
- **Information about fire - and explosion protection:** No special measures required.

(Contd. on page 3)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 2)

- 7.2 Conditions for safe storage, including any incompatibilities
- Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- 7.3 Specific end use(s) : No further relevant information available.

**SECTION 8: Exposure controls/personal protection**

- Additional information about design of technical facilities: No further data; see section 7.

· 8.1 Control parameters

- Ingredients with limit values that require monitoring at the workplace:

**CAS: 14808-60-7 Quartz (SiO<sub>2</sub>) (25-50%)**

MAK (Switzerland)	Long-term value: 0.15 a mg/m <sup>3</sup> P C1a SSc;
-------------------	---

**CAS: 56-81-5 glycerol (10-25%)**

WEL (Great Britain)	Long-term value: 10 mg/m <sup>3</sup>
MAK (Switzerland)	Short-term value: 100 e mg/m <sup>3</sup> Long-term value: 50 e mg/m <sup>3</sup> SSc;

**CAS: 64-17-5 ethanol (≤2.5%)**

WEL (Great Britain)	Long-term value: 1920 mg/m <sup>3</sup> , 1000 ppm
MAK (Switzerland)	Short-term value: 1920 mg/m <sup>3</sup> , 1000 ppm Long-term value: 960 mg/m <sup>3</sup> , 500 ppm SSc;

**CAS: 9005-25-8 Starch (<1%)**

WEL (Great Britain)	Long-term value: 10* 4** mg/m <sup>3</sup> *total inhalable **respirable
MAK (Switzerland)	Long-term value: 3 a mg/m <sup>3</sup>

**CAS: 1310-58-3 potassium hydroxide (<0.1%)**

WEL (Great Britain)	Short-term value: 2 mg/m <sup>3</sup>
MAK (Switzerland)	Long-term value: 2 e mg/m <sup>3</sup>

· DNELs

**CAS: 64-17-5 ethanol**

Dermal	DNEL Worker long-term systemic	343 mg/kg bw/d
Inhalative	DNEL Worker long-term systemic	950 mg/m <sup>3</sup>
	DNEL Worker acute/short-term local effects	19 mg/m <sup>3</sup>

· PNECs

**CAS: 64-17-5 ethanol**

Oral	PNEC sediment (fresh water)	3.6 mg/kg
	PNEC STP	0.72 mg/L
	PNEC aqua (fresh water)	0.96 mg/L
	PNEC aqua (marine water)	0.79 mg/L
	PNEC sediment (marine water)	2.9 mg/kg

(Contd. on page 4)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 3)

PNEC soil	0.63 mg/kg
-----------	------------

- **Additional information:**  
The most current valid lists have been used as a basis for the production of this document.
- **8.2 Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required.
- **Protection of hands:** Not required.
- **Gloves material :** -
- **Eye protection:** Not required.

**SECTION 9: Physical and chemical properties**

- **9.1 Information on basic physical and chemical properties**
- **General Information**
- **Appearance:**
  - **Form:** Pasty
  - **Colour:** Pink
  - **Odour:** Characteristic
  - **Odour threshold:** Not determined.
- **pH-value:** Not applicable.
- **Change in condition**
  - **Melting point/freezing point:** Undetermined.
  - **Initial boiling point and boiling range:** Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gas):** Not determined.
- **Ignition temperature:** 400 °C
- **Decomposition temperature:** Not determined.
- **Auto-ignition temperature:** Product is not self igniting.
- **Explosive properties:** Product does not present an explosion hazard.
- **Explosion limits:**
  - **Lower:** 0,9 Vol %
  - **Upper:** Not determined.
- **Vapour pressure at 20 °C:** 0,1 hPa
- **Density:** Not determined.
- **Relative density :** Not determined.
- **Vapour density :** Not applicable.
- **Evaporation rate :** Not applicable.
- **Solubility in / Miscibility with water:** Insoluble.
- **Segregation coefficient (n-octanol/water):** Not determined.
- **Viscosity:**
  - **Dynamic:** Not applicable.

(Contd. on page 5)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 4)

**Kinematic:** Not applicable.

· **Solvent content:**

**VOC (EC)** 0.00 %

· **9.2 Other information :** No further relevant information available.

**SECTION 10: Stability and reactivity**

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

**SECTION 11: Toxicological information**

- **11.1 Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

**CAS: 10101-52-7 Zirconium(IV) silicate**

Inhalative	LC50/4 h	1.5 mg/l (ATE)
------------	----------	----------------

**CAS: 64-17-5 ethanol**

Oral	LD50	7,060 mg/kg (rat)
------	------	-------------------

Inhalative	LC50/4 h	20,000 mg/l (rat)
------------	----------	-------------------

- **Primary irritant effect:**
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

**SECTION 12: Ecological information**

· **12.1 Toxicity**

· **Aquatic toxicity:**

**CAS: 64-17-5 ethanol**

EC50 72h	275 mg/l (Algae)
----------	------------------

EC50 48h	12,340 mg/l (daphnia)
----------	-----------------------

EC50 96h	13,000 mg/l (fish)
----------	--------------------

· **12.2 Persistence and degradability :** No further relevant information available.

(Contd. on page 6)



**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 5)

- **12.3 Bioaccumulative potential** : No further relevant information available.
- **12.4 Mobility in soil** : No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **12.5 Results of PBT and vPvB assessment**
- **PBT**: Not applicable.
- **vPvB**: Not applicable.
- **12.6 Other adverse effects** : No further relevant information available.

**SECTION 13: Disposal considerations**

- **13.1 Waste treatment methods**
- **Recommendation** : Smaller quantities can be disposed of with household waste.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

**SECTION 14: Transport information**

- |   |   |
|---|---|
| · <b>14.1 UN-Number</b><br>· <b>ADR, ADN, IMDG, IATA</b>                                    | Not classified as hazardous for transport |
| · <b>14.2 UN proper shipping name</b><br>· <b>ADR, ADN, IMDG, IATA</b>                      | Not classified as hazardous for transport |
| · <b>14.3 Transport hazard class(es)</b><br>· <b>ADR, ADN, IMDG, IATA</b><br>· <b>Class</b> | Not classified as hazardous for transport |
| · <b>14.4 Packing group</b><br>· <b>ADR, IMDG, IATA</b>                                     | Not classified as hazardous for transport |
| · <b>14.5 Environmental hazards:</b><br>· <b>Marine pollutant:</b>                          | No  |
| · <b>14.6 Special precautions for user</b>  | Not applicable.                           |
| · <b>14.7 Transport in bulk according to Annex II of Marpol and the IBC Code</b>            | Not applicable.                           |
| · <b>UN "Model Regulation":</b>   | Not classified as hazardous for transport |

**SECTION 15: Regulatory information**

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.

(Contd. on page 7)

GB

**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: DETARTRINE Z**

(Contd. of page 6)

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

**SECTION 16: Other information**

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

**· Relevant phrases**

*H225 Highly flammable liquid and vapour.*

*H315 Causes skin irritation.*

*H319 Causes serious eye irritation.*

*H332 Harmful if inhaled.*

*H335 May cause respiratory irritation.*

**· Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*IATA: International Air Transport Association*

*GHS: Globally Harmonised System of Classification and Labelling of Chemicals*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*VOC: Volatile Organic Compounds (USA, EU)*

*DNEL: Derived No-Effect Level (REACH)*

*PNEC: Predicted No-Effect Concentration (REACH)*

*LC50: Lethal concentration, 50 percent*

*LD50: Lethal dose, 50 percent*

*PBT: Persistent, Bioaccumulative and Toxic*

*SVHC: Substances of Very High Concern*

*vPvB: very Persistent and very Bioaccumulative*

*Flam. Liq. 2: Flammable liquids – Category 2*

*Acute Tox. 4: Acute toxicity – Category 4*

*Skin Irrit. 2: Skin corrosion/irritation – Category 2*

*Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

*STOT SE 3: Specific target organ toxicity (single exposure) – Category 3*