

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

- **1.1 Product identifier**
- **Trade name:** ENDOSOLV
- **Article number / formula:** 279S
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the preparation :** dental product
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
SEPTODONT.  
58 rue du Pont de Creteil  
94100 Saint Maur des Fosses  
France  
Téléphone: +33(0)149767000  
Fax: +33(0)148855401  
msds@septodont.com  
http://www.septodont.com
- **1.4 Emergency telephone number:**  
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
Belgique - Tel : 32 070/245 245  
UK : +44 (0)-870-190-6777 (CARECHEM24)

**SECTION 2: Hazards identification**

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**



Flam. Liq. 2 H225 Highly flammable liquid and vapour.



Eye Irrit. 2 H319 Causes serious eye irritation.

STOT SE 3 H336 May cause drowsiness or dizziness.

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008**  
The product is classified and labelled according to the CLP regulation.
- **Hazard pictograms**



GHS02 GHS07

- **Signal word** Danger
- **Hazard-determining components of labelling:**  
ethyl acetate

(Contd. on page 2)

GB

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 1)

**· Hazard statements**

- H225 Highly flammable liquid and vapour.
- H319 Causes serious eye irritation.
- H336 May cause drowsiness or dizziness.

**· Precautionary statements**

- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P280 Wear eye protection / face protection.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P403+P233 Store in a well-ventilated place. Keep container tightly closed.
- P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

**· Labelling of packages where the contents do not exceed 125 ml**

**· Hazard pictograms**



GHS02 GHS07

**· Signal word Danger**

**· Hazard-determining components of labelling:**

ethyl acetate

**· Hazard statements** Not classified as hazardous

**· Precautionary statements**

- P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
- P280 Wear eye protection / face protection.
- P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P403+P233 Store in a well-ventilated place. Keep container tightly closed.

**· 2.3 Other hazards**

**· Results of PBT and vPvB assessment**

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

**SECTION 3: Composition/information on ingredients**

**· 3.2 Chemical characterisation: Mixtures**

· **Description:** Mixture: consisting of the following components.

**· Dangerous components:**

CAS: 141-78-6 EINECS: 205-500-4 Index number: 607-022-00-5 RTECS: AH 5425000 Reg.nr.: 1-2119475103-46-XXXX	ethyl acetate Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336	50-100%
CAS: 84145-37-9 EINECS: 282-263-3 Index number: 607-130-00-2	2(or 3)-methylbutyl acetate Flam. Liq. 3, H226	5-10%

**· SVHC** No

(Contd. on page 3)

## Safety data sheet

According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 2)

- **Additional information:** For the wording of the listed risk phrases, refer to section 16.

### SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **After excessive inhalation:** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult a doctor.
- **4.2 Most important symptoms and effects, both acute and delayed :**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed :**  
No further relevant information available.

### SECTION 5: Firefighting measures

- **5.1 Suitable extinguishing agents:**  
CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

### SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :** Not required.
- **6.2 Environmental precautions:** Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).  
Dispose contaminated material as waste according to item 13.
- **6.4 Reference to other sections :**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protective equipment.  
See Section 13 for disposal information.

### SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :** No special precautions are necessary if used correctly.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **7.3 Specific end use(s) :** No further relevant information available.

### SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see section 7.

(Contd. on page 4)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 3)

· **8.1 Control parameters**

· **Ingredients with limit values that require monitoring at the workplace:**

**CAS: 141-78-6 ethyl acetate (50-100%)**

WEL (Great Britain)	Short-term value: 400 ppm Long-term value: 200 ppm
IOELV (EU)	Short-term value: 1468 mg/m <sup>3</sup> , 400 ppm Long-term value: 734 mg/m <sup>3</sup> , 200 ppm
MAK (Switzerland)	Short-term value: 2800 mg/m <sup>3</sup> , 800 ppm Long-term value: 1400 mg/m <sup>3</sup> , 400 ppm SSc;

· **DNELs**

**CAS: 141-78-6 ethyl acetate**

Oral	DNEL Public long-term systemic	4.5 mg/kg bw/d
Dermal	DNEL Public long-term systemic	37 mg/kg bw/d
	DNEL Worker long-term systemic	63 mg/kg bw/d
Inhalative	DNEL Worker long-term local effects	734 mg/m <sup>3</sup>
	DNEL Public acute/short-term local effects	734 mg/m <sup>3</sup>
	DNEL Worker acute/short-term local effects	1,468 mg/m <sup>3</sup>

· **PNECs**

**CAS: 141-78-6 ethyl acetate**

PNEC STP	mg/L (650)
PNEC aqua (fresh water)	1.125 mg/L
PNEC aqua (intermittent releases)	1.65 mg/L
PNEC aqua (marine water)	0.026 mg/L
PNEC sediment (marine water)	0.125 mg/kg
PNEC soil	0.24 mg/kg

· **Additional information:**

The most current valid lists have been used as a basis for the production of this document.

· **8.2 Exposure controls**

· **Personal protective equipment:**

· **General protective and hygienic measures:**

Keep away from food, beverages and petfood.  
Wash hands before breaks and at the end of work.

· **Respiratory protection:** Not required.

· **Protection of hands:**



Protective gloves

· **Gloves material :** Butyl rubber, BR

(Contd. on page 5)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 4)

· **Eye protection:**



Tightly sealed goggles

**SECTION 9: Physical and chemical properties**

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· <b>Form:</b>	Liquid
· <b>Colour:</b>	Colourless
· <b>Odour:</b>	Characteristic
· <b>Odour threshold:</b>	Not determined.

· **pH-value:** Not determined.

· **Change in condition**

· <b>Melting point/freezing point:</b>	Undetermined.
· <b>Initial boiling point and boiling range:</b>	Undetermined.

· **Flash point:** -1 °C

· **Flammability (solid, gas):** Not applicable.

· **Ignition temperature:** Not applicable.

· **Decomposition temperature:** Not determined.

· **Auto-ignition temperature:** Product is not self igniting.

· **Explosive properties:** Product does not present an explosion hazard.

· **Explosion limits:**

· <b>Lower:</b>	Not determined.
· <b>Upper:</b>	Not determined.

· **Vapour pressure:** Not determined.

· **Density:** Not determined.

· **Relative density :** Not determined.

· **Vapour density :** Not determined.

· **Evaporation rate :** Not determined.

· **Solubility in / Miscibility with water:**

Not miscible or difficult to mix.

· **Segregation coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

· **Dynamic:** Not determined.

· **Kinematic:** Not determined.

· **9.2 Other information :** No further relevant information available.

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 5)

**SECTION 10: Stability and reactivity**

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

**SECTION 11: Toxicological information**

- **11.1 Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

**CAS: 141-78-6 ethyl acetate**

Oral	LD50	5,620 mg/kg (rabbit)
Inhalative	LC50/4 h	1,600 mg/l (rat)

- **Primary irritant effect:**
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:**  
Causes serious eye irritation.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure**  
May cause drowsiness or dizziness.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

**SECTION 12: Ecological information**

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability :** No further relevant information available.
- **12.3 Bioaccumulative potential :** No further relevant information available.
- **12.4 Mobility in soil :** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects :** No further relevant information available.

GB

(Contd. on page 7)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 6)

**SECTION 13: Disposal considerations**

- 13.1 Waste treatment methods
- Recommendation :



Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- European waste catalogue

HP 3	Flammable
HP 4	Irritant - skin irritation and eye damage
HP 5	Specific Target Organ Toxicity (STOT)/Aspiration Toxicity

- Uncleaned packaging:
- Recommendation: Disposal must be made according to official regulations.

**SECTION 14: Transport information**

- 14.1 UN-Number
- ADR, IMDG, IATA

UN1993

- 14.2 UN proper shipping name
- ADR

UN1993 FLAMMABLE LIQUID, N.O.S. (ETHYL ACETATE), special provision 640D  
FLAMMABLE LIQUID, N.O.S. (ETHYL ACETATE)

- IMDG, IATA

- 14.3 Transport hazard class(es)

- ADR



- Class
- Label

3 (F1) Flammable liquids.  
3

- IMDG, IATA



- Class
- Label

3 Flammable liquids.  
3

- 14.4 Packing group
- ADR, IMDG, IATA

II

- 14.5 Environmental hazards:
- Marine pollutant:

Yes

- 14.6 Special precautions for user
- Danger code (Kemler):

Warning: Flammable liquids.  
33

(Contd. on page 8)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 7)

· <b>EMS Number:</b>	<i>F-E,S-E</i>
· <b>Stowage Category</b>	<i>B</i>
· <b>14.7 Transport in bulk according to Annex II of Marpol and the IBC Code</b>	<i>Not applicable.</i>
· <b>Transport/Additional information:</b>	
· <b>ADR</b>	
· <b>Limited quantities (LQ)</b>	<i>1L</i>
· <b>Excepted quantities (EQ)</b>	<i>Code: E2</i> <i>Maximum net quantity per inner packaging: 30 ml</i> <i>Maximum net quantity per outer packaging: 500 ml</i>
· <b>Transport category</b>	<i>2</i>
· <b>Tunnel restriction code</b>	<i>D/E</i>
· <b>IMDG</b>	
· <b>Limited quantities (LQ)</b>	<i>1L</i>
· <b>Excepted quantities (EQ)</b>	<i>Code: E2</i> <i>Maximum net quantity per inner packaging: 30 ml</i> <i>Maximum net quantity per outer packaging: 500 ml</i>
· <b>UN "Model Regulation":</b>	<i>UN 1993 FLAMMABLE LIQUID, N.O.S., SPECIAL PROVISION 640D (ETHYL ACETATE), 3, II</i>

**SECTION 15: Regulatory information**

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **Seveso category P5c** FLAMMABLE LIQUIDS
- **Qualifying quantity (tonnes) for the application of lower-tier requirements** 5,000 t
- **Qualifying quantity (tonnes) for the application of upper-tier requirements** 50,000 t
- **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

**SECTION 16: Other information**

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

- **Relevant phrases**
- H225 Highly flammable liquid and vapour.*
- H226 Flammable liquid and vapour.*
- H319 Causes serious eye irritation.*
- H336 May cause drowsiness or dizziness.*

· **Classification according to Regulation (EC) No 1272/2008**

<i>Flammable liquids</i>	<i>Bridging principles</i>
<i>Serious eye damage/eye irritation</i> <i>Specific target organ toxicity (single exposure)</i>	<i>The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.</i>

(Contd. on page 9)





**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: ENDOSOLV**

(Contd. of page 8)

**· Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*

*IMDG: International Maritime Code for Dangerous Goods*

*IATA: International Air Transport Association*

*GHS: Globally Harmonised System of Classification and Labelling of Chemicals*

*EINECS: European Inventory of Existing Commercial Chemical Substances*

*ELINCS: European List of Notified Chemical Substances*

*CAS: Chemical Abstracts Service (division of the American Chemical Society)*

*DNEL: Derived No-Effect Level (REACH)*

*PNEC: Predicted No-Effect Concentration (REACH)*

*LC50: Lethal concentration, 50 percent*

*LD50: Lethal dose, 50 percent*

*PBT: Persistent, Bioaccumulative and Toxic*

*SVHC: Substances of Very High Concern*

*vPvB: very Persistent and very Bioaccumulative*

*Flam. Liq. 2: Flammable liquids – Category 2*

*Flam. Liq. 3: Flammable liquids – Category 3*

*Eye Irrit. 2: Serious eye damage/eye irritation – Category 2*

*STOT SE 3: Specific target organ toxicity (single exposure) – Category 3*

GB