

NEW Septodont
**hydraulic temporary
restorative
material**



Cimpat W
Cimpat P



Cimpat W & Cimpat P

Hydraulic temporary restorative material


Intended to be used as a temporary filling material. Hydraulic temporary restorative materials are based on calcium sulfate: they set upon contact with the saliva in the oral cavity⁽¹⁾

Composition

Co-polymer, Zinc oxide, Calcium sulfate

Features	Advantages	Benefits
Paste texture	Easy application and removal	Dentist time optimization
Hydraulic temporary restorative material	Hermetic seal of the cavity Excellent stability when set	Effective and easy temporary filling
Two colors available	Can be used for anterior and posterior teeth	Convenience to use

- Sets upon the contact with oral fluid
- Setting time: 30 min⁽²⁾



(1) Thu, Khaing & Aye, Khin & Htang, Aung. (2013). In vitro study of Coronal Leakage of Four Temporary Filling Material

(2) Test certificate #16-02-005 issued by Dental Devices Testing & Evaluation Center Yonsei University, College of Dentistry

Distributor:

Septodont - 58 rue du Pont de Créteil - 94107 Saint-Maur-des-Fossés Cedex - France

Tél. : +33 (0)1 49 76 70 00 - Fax : +33 (0)1 48 85 54 01

Please visit our website for more information: www.septodont.com

